

Retainers


AOA employs various methods for fabricating **Hawley** retainers, with the most favored approach being the submission of digital scans prior to the debonding of brackets. This process allows AOA ample time for fabrication, ensuring that patients receive their retainer on debonding day. Many practitioners have found this method to be an efficient means of transitioning from active orthodontic treatment to lasting retention. To utilize this method, it is necessary to take scans directly over the brackets. Additionally, to meet our customers' needs, AOA still accepts plaster or stone models for retainer fabrication.

AOA's retainer department can process virtually any clasping request; finger springs, and acrylic modifications such as bite planes and occlusal coverage. A wide variety of expansion and single tooth movement screws are available.


Anterior and posterior pontics may be requested in most of our designs, or in your custom design. Just indicate shade preference and we can match the shade.. Tissue compatible flesh or clear opaque acrylic may be requested for blending of the pontic against the patient's contour.



Invisible



QCM



Wrap-around Hawley



Custom Acrylics



Hawley

AOA

001-805 Rev B